

## MONTGOMERY® ESOPHAGEAL TUBE

### DESCRIPTION:

The **Montgomery® Esophageal Tube** is used to bridge the gap between the pharyngostome and esophagostome following laryngoesophagectomy and first stage reconstruction of the cervical esophagus. The tube features a funnel-shaped superior end to conform to the hypopharynx, and two firm oval bands designed to ensure proper positioning and to prevent leakage of saliva after laryngoesophagectomy.

The Montgomery® Esophageal Tube has three functions: 1) it tends to eliminate the problem of profuse salivary leakage by way of the pharyngostome, 2) it maintains a widely patent pharyngostome and esophagostome, and, 3) it creates a trough between the pharyngostome and esophagostome facilitating the second stage procedure.

### INDICATIONS FOR USE:


For two-stage reconstruction of the cervical esophagus.

### SUPPLIED:

Sterile, 1 tube per carton

### PRODUCT CODES:

**322000**    **ESOPHAGEAL TUBE, CLEAR**  
**322000R**    **ESOPHAGEAL TUBE, RADIOPAQUE**



### FEATURES

- Implant grade silicone.
- Non-adherent smooth surface.
- Available in clear or radiopaque material.

