

OTOSILK™ GRAFT DRESSING

DESCRIPTION:

The OTOSILK™ Graft Dressing is a 4" x 6" sheet of blue, non-adherent material for use as a carrier for split thickness skin grafts placed in the external auditory canal. The graft is placed on the dressing and trimmed to the desired size and shape. The dressing and graft combination is then placed in the canal. Following the healing process, the graft dressing is removed, leaving the newly healed graft in position.

INDICATIONS FOR USE:

For use as a temporary dressing in the external auditory canal following surgery.

SUPPLIED:

Sterile, 10 dressings per carton
(1 dressing / sterile pouch; 10 pouches / dispenser carton)

PRODUCT CODE:

OT-4601 OTOSILK™ GRAFT DRESSING



FEATURES

- 4" x 6" sheet, easily trimmed.
- Non-adherent material.
- Blue color for increased visualization.