

## OTOSILK™ OTOLOGICAL STRIPS

### DESCRIPTION:

OTOSILK™ **Otological Strips** are used in the process of packing the external auditory canal following surgical procedures such as tympanoplasty or stapedectomy. The strips can be moistened prior to use using either sterile saline or medication solution. The blue color of the strips aids in visualization during removal. The ends of the OTOSILK™ strips have been heat-sealed to prevent unraveling or shedding. Each strip measures .25" wide by 1.75" long.

### INDICATIONS FOR USE:


For use as a temporary dressing in the external auditory canal following surgery.

### SUPPLIED:

Sterile, 200 strips per carton  
(10 strips / pouch, 20 pouches / dispenser carton)

### PRODUCT CODE:

**OT-17510 OTOSILK™ OTOLOGICAL STRIPS**



### FEATURES

- Convenient, economical pre-cut strips.
- Non-adherent material.
- Blue color for increased visualization.
- Heat-sealed edges.