

## SALMAN FES STENT

### DESCRIPTION:

The **Salman FES Stent** is designed to help prevent potential complications resulting from functional endoscopic sinus surgery (FES). The stent has an oval-shaped base that may be trimmed for different anatomies and is designed to fit between the middle turbinate and lateral nasal wall. The two “fingers” arising from the base are introduced into the maxillary sinus through the middle meatal antrostomy. These “fingers” will prevent dislodgment of the stent, thus eliminating the need for any suturing.

### INDICATIONS FOR USE:

Prevention of post-operative adhesions.  
Prevention of narrowing or closure of the middle meatal antrostomy.

### SUPPLIED:

Sterile, 1 stent per carton

### PRODUCT CODE:

**SS2000**     **SALMAN FES STENT**



### FEATURES

- Medical grade silicone material.
- White color for maximum visualization.
- Smooth non-adherent surface.
- Can be trimmed prior to insertion.

