

## EPI-MAX™ EPISTAXIS CATHETER

### DESCRIPTION:

The **Epi-Max™** is a two-balloon epistaxis catheter designed to control intranasal bleeding. As a nasal packing, it can be used following septoplasty, rhinoplasty, and other intra-nasal surgical procedures. For epistaxis, the smaller balloon controls posterior bleeding, while the larger balloon handles anterior bleeding. The included syringe is used to expand the balloons with normal saline. The volume of saline injected determines the ultimate size and pressure of the balloons. The result is a controlled uniform pressure to manage bleeding. The integral airway tube provides a passage for breathing while the balloons are inflated.

### INDICATIONS FOR USE:

For control of anterior and/or posterior nasal epistaxis.  
For nasal packing following septoplasty, rhinoplasty, and other intra-nasal surgical procedures.

### SUPPLIED:

Sterile, 1 catheter per carton

### PRODUCT CODE:

**60-402 EPI-MAX™ EPISTAXIS CATHETER**



### FEATURES

- Medical grade silicone with blue tint.
- Dual balloons with adjustable size and pressure.
- Integral airway.
- 20cc syringe included.