

HAR-EL PHARYNGEAL TUBE

DESCRIPTION:

The **Har-El Pharyngeal Salivary Bypass Tube** allows the surgeon to handle a wide variety of indications for treatment of esophageal fistulae. The tube features a unique superior funnel shape designed for anchoring at the level of the tongue base. The tube has applications for: 1) the treatment of orocutaneous and pharyngocutaneous fistulae following head and neck oncologic resection, 2) placement in patients with an existing fistula prior to reconstructive procedures, and, 3) prophylactic placement in patients who are at high risk for development of a fistula.

INDICATIONS FOR USE:

Treatment of esophageal fistulae.

SUPPLIED:

Sterile, 1 tube per carton

PRODUCT CODES:

355120 **SIZE 12mm, CLEAR**
355120R **SIZE 12mm, RADIOPAQUE**



FEATURES

- Implant grade silicone.
- Non-adherent smooth surface.
- Available in clear or radiopaque material.

